



### Suggested Printers and DVRs

Model	Type	Dimension W,H,D, [mm]	Weight [Kg]	Interfaces
HP Color LaserJet Pro M454dw	Laser color	412, 295, 469	18.9	LAN, WI-FI
Canon Pixma TR150	Inkjet	322, 66, 185	2.1	USB, WI-FI
Epson Workforce WF-110W	Inkjet	309, 61, 159	1.7	USB, WI-FI
Sony UP-D711MD	Thermal B/W	140, 70, 125	1.4	USB
Sony UP-X898MD	Thermal B/W	154, 88, 240	2.5	Composite, USB
Sony UP-D25MD	Thermal color	212, 398, 98	5.7	USB
Sony UP-DR80MD	A4 Thermal color	317, 207, 425	7.3	USB
Sony HVO-550MD	DVD Recorder	212, 105, 288	3.2	S-Video, HDMI

Please note that the Ultrasound Peripherals are system/configuration dependent.

Please check with your Service representative the compatibility between your MyLab™ system and the desired printer or DVR, and ask for the installation instructions and certified drivers to use.

Installing with wrong drivers could compromise the correct system functionality.

Technology and features are system/configuration dependent. Specifications subject to change without notice.

Information might refer to products or modalities not yet approved in all countries.

For further details, please contact your Esaote sales representative: [info\(at\)esaote.com](mailto:info(at)esaote.com)