



### Suggested Printers and DVRs

Model	Type	Dimension W,H,D, [mm]	Weight [Kg]	Interfaces
HP Color LaserJet Pro 4202dw	Laser color	421, 264, 427	16.4	LAN, WI-FI
Canon Maxify BX110	Inkjet	322, 66, 210	2.3	WI-FI
Canon Pixma TR160	Inkjet	322, 66, 210	2.3	WI-FI
Sony UP-D711MD	Thermal B/W	140, 70, 125	1.4	USB
Sony UP-D899MD	Thermal B/W	154, 88, 208	2.0	USB
Sony UP-D25MD	Thermal color	212, 398, 98	5.7	USB
Sony UP-DR80MD	A4 Thermal color	317, 207, 425	7.3	USB
Sony HVO-550MD	DVD Recorder	212, 105, 288	3.2	S-Video, HDMI

Please note that the Ultrasound Peripherals are system/configuration dependent.

Please check with your Service representative the compatibility between your MyLab™ system and the desired printer or DVR, and ask for the installation instructions and certified drivers to use. Installing with wrong drivers could compromise the correct system functionality.

Technology and features are system/configuration dependent. Specifications subject to change without notice.  
Information might refer to products or modalities not yet approved in all countries.  
For further details, please contact your Esaote sales representative: [info@esaote.com](mailto:info@esaote.com)